

Adscope™ 608 Convertible Clinician Stethoscope

ADSCOPE™



AFD multifrequency diaphragm

Convert pediatric diaphragm to adult bell in seconds

Accessory kit with extra Adsoft™ Plus eartips

Combining AFD technology with a convertible design makes the new Adscope™ 608 the most versatile of our Clinician scopes.

The Adscope™ 608 series features:

- Convertible chestpiece provides adult diaphragm/bell or adult/pediatric diaphragm configurations
- Combination diaphragm/bell adult chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability
- ADC's AFD technology provides enhanced acoustic response: use light pressure for low, bell-type frequency, firm pressure for high, diaphragm-type frequency
- Extra-large bell (1 3/8" diameter) for unsurpassed low-frequency response
- Ultra-sensitive diaphragm (1 3/4" diameter) for greater amplification and crisper high-frequency response
- Non-chill bell and diaphragm retaining rims for patient comfort
- Clinician headset boasts reinforcing yoke molded into flexible 22" PVC tubing and stainless steel binaurals fixed at 15° angle for greater comfort
- Adsoft™ Plus silicone eartips with new snap-on connections for the ultimate in wearing comfort and acoustic seal
- Includes accessory kit with two pairs of spare Adsoft™ Plus eartips
- Scope ID tag included
- Weighs 5.77 oz. (165g)
- Overall length 31"
- Assembled, inspected, and packaged in the USA
- Lifetime warranty covers all parts

Color Choices for 608**

BK	BD	DG	FG	FL	FRB	FS	FV
G	LV	MCB	MCA	MG	MRS	N	O
P	PBCA	R	RB	TL	TQ	V	

NOTE: To ensure correct color when ordering. Add color suffix to model number (e.g., 608DG)

Ordering Information

Item	Description
608**	Adscope™ 608 Stainless Steel, Adult, (Specify Color)

